



For guidance in other languages and a video guide visit: www.bit.ly/fitlondon or scan this QR code

PATIENT INSTRUCTIONS FOR COLLECTING YOUR SAMPLE

BEFORE YOU BEGIN

You have been asked to complete this test as soon as possible.


Use a black or blue pen to write **YOUR NAME**, **YOUR DATE OF BIRTH** and **NHS HOSPITAL NUMBER** on the label.

FIT (Faecal Immunochemical Test) is an investigation that looks for the presence of blood in the stools (poo) when it is too small an amount to be visible.

Your clinician has ordered this test for you as part of the investigation of your symptoms to rule out cancer.


The FIT kit includes a small plastic container. You should complete the test and return this to your GP practice or drop off location below as soon as possible. Please follow the step-by-step instructions below to collect your stool (poo).

1. Label your sample



Write your **FIRST NAME** and **SURNAME (N)**, **DATE OF BIRTH** not age (A) and **NHS HOSPITAL NUMBER** along with the **DATE (D)** you collected your stool (poo) where indicated using a black or blue pen. If the label already has your details on it please check they are correct.

3. Take a sample of your stool (poo) with the extraction device



Unscrew the top of the sample tube anti-clockwise, remove the stick and **scrape the end of the stick along the stool (poo)**. Only a very small amount of the stool (poo) is needed on the stick.

2. Ideas about how to catch the stool (poo) sample




It is important that your stool (poo) sample does not touch the toilet water.

There are different ways to collect your poo, try:

- Hold folded lengths of toilet paper in your hand.
- Put a small, disposable plastic bag/glove over your hand.
- Poo into a clean, disposable container lined with a sheet of toilet paper so that you can easily slip the stool (poo) into the toilet after.

You may find a way that is easier for you, but make sure you catch your stool (poo) sample before it touches the water.

4. Close and seal the container



Replace the sample stick back into the tube. Make sure that the tube clicks closed and place it in the plastic bag with your request form.

Don't forget to wash your hands!

6. Results

The hospital will process the sample and communicate the results to the appropriate clinician or your GP and yourself within 10 days. Make a note on your calendar, if you have not heard from us in two weeks, please contact FiT@Bhrut.nhs.uk or call **01708 435 036**. If blood is found in your poo this does not necessarily mean you have cancer, but you will need to be seen for further investigations. For further information about bowel cancer please contact:


Bowel Cancer UK:

<https://www.bowelcanceruk.org.uk/about-bowel-cancer/diagnosis/visiting-your-gp/>

Macmillan Cancer Support 0808 808 0000

Cancer Research UK 0808 800 4040

5. Hand in your sample



Take your sample to your GP practice or drop off location below. Please submit the sample **within 7 days**.



WHAT TO DO WITH YOUR SAMPLE

Please drop your sample at **your GP surgery** or at one of the following locations:

Queen's Hospital: Outpatients Team 3 & 4

King George Hospital: Pathology Unit (Cedar Centre)

Havering: Cranham Health Centre (RM14 1RG), Elm Park Health Centre (RM12 4YG), Raphael House/Victoria Centre (RM1 4HP), Harold Hill Health Centre (RM3 9SU)

Barking and Dagenham: Chadwell Heath Health Centre (RM6 6RT), Porters Avenue Health Centre (RM8 2EQ), Barking Hospital (IG11 9LX) and Thames View Health Centre (IG11 0LG)

Redbridge: Hainault Health Centre (IG7 4DF), Loxford Polyclinic (IG1 2SN), Barley Court Clinic (IG3 8XJ).